


15 & 50mL Polypropylene Storage Boxes

PolarSafe polypropylene freezer boxes can accommodate 15&50mL centrifuge tubes for long term storage; an excellent alternative to cardboard storage boxes. Boxes can be stored to -90°C (130°F) and are autoclavable to 122°C (251.6°F). (Dimensions: 5 3/4L x 5 3/4W x 4 7/8"H)



- Clear lids for easy visual inspection
- Easy-to-read numerical grid for tube reference and identification
- White top portion allows for easy labeling or marking
- Accommodates most 4" stainless steel freezer racks

Catalog No.	Description	Qty.
R0015	36-Place PP Freezer Box for 15mL Tubes	2/Pack
R0050	16-Place PP Freezer Box for 50mL Tubes	2/Pack

15 & 50mL Cardboard Storage Box

Cardboard freezer boxes are constructed of sturdy, durable, white cardboard with a moisture repellant coating that protects the box from the affects of liquid nitrogen or condensation from refrigerators or freezers. Boxes feature small drain holes on the bottom for safe liquid nitrogen drainage or to allow excess moisture to escape.

Boxes and dividers sold separately.



Tubes sold separately


- Moisture repellant coating prevents water or liquid nitrogen from saturating box
- Easy to write on surface for identification
- Drain holes on bottom - allows for safe liquid nitrogen drainage and prevents excess moisture build-up

Catalog No.	Description	Qty.
R3014A	Cardboard Freezer Box (5 3/4 x 5 3/4 x 4 7/8 in.)	1/Each
R3022	36-Place Divider for 15mL Tubes	12/Pack
R3021	16-Place Divider for 50mL Tubes	12/Pack



15 & 50mL Polypropylene Tubes

Disposable polypropylene tubes and caps with recessed lid and sealing ring to prevent leaking. Frosted writing surface allows for permanent identification and labeling. Autoclavable to 121°C and can be stored in freezers to -80°C. Manufactured in ISO 9001:2000 and 100,000 grade cleanrooms. Sterile.



- Easy-to-read black graduations are accurate to $\pm 2\%$, 1mL increments (15mL) or 2.5mL increments (50mL)


- Available in cardboard racks or in resealable zip closure bags

- Can be spun up to 12,000g (15mL) & 15,000g (50mL)

- Conical Bottom

- Sterilized with gamma irradiation, RNase & DNase Free, DNA free, Non-pyrogenic

- Chloroform resistant graduations and writing area



Environmentally friendly, fully recyclable cardboard racks

Catalog No.	Capacity	Packaging	Qty.
TC1500	15mL	Resealable Bag	25/Bag, 500/Case
TC1500R	15mL	Cardboard Rack	25/Rack, 500/Case
TC5000	50mL	Resealable Bag	25/Bag, 500/Case
TC5000R	50mL	Cardboard Rack	25/Rack, 500/Case



Ph: 800-886-8675 • Fx: 847-622-0454
www.argos-tech.com

Distributed by: