




# Argos *BackSpring*

## 1.5mL Microcentrifuge Tubes


Argos BackSpring tubes are designed to reduce repetitive stress syndrome caused by repeatedly opening and closing microcentrifuge tubes. Unique design allows the tube to open from the back, versus the traditional front open style.



**Traditional front open style tubes can leave thumbs sore and stressed after repeated opening.**



**The BackSpring's™ unique design allows the tube to snap open when pressure is applied to the back of cap.**



**Enables quick, comfortable opening and closing of lid; eliminating repetitive stress to thumbs.**

- Manufactured of homopolymer plastic
- RNase / DNase Free
- Withstands freezing to -90°C
- Boil-Proof
- Easy to read graduations / frosted writing surface
- Can be spun up to 30,000g

#### Ordering Info:

Cat. No.	Description	Qty.
T0015	1.5mL BackSpring Tubes	500/Pk.

**ARGOS**  
TECHNOLOGIES

Ph: 800-886-8675 • Fx: 847-622-0454  
www.argos-tech.com

Distributed by: |