



DR NAsE FREE



- *Ideal for cleaning labware that cannot be autoclaved*
- *Residue Free*
- *Safe solution for eliminating RNase and DNase from laboratory surfaces*

DR NAsE Free is a surface decontaminant that is specially formulated to eliminate RNase and DNase from laboratory glassware, plastics, stainless steel and any laboratory surface. Ideal for cleaning equipment, pipettors, gel boxes and any labware that cannot be autoclaved. Simply apply DR NAsE Free to area to be decontaminated and rinse. DR NAsE Free leaves no residue that could interfere with subsequent enzymatic reactions. DR NAsE Free offers a non-toxic formula that replaces carcinogenic chemicals such as DEPC. No need to bake glassware. The product is non-abrasive and leaves no residue that can interfere with gel polymerization. The product is non-biohazardous and can be stored on the laboratory bench.

Cat. No.	Size Specifications	Qty.	List Price
D6002	720 mL / spray Bottle	CS/12	\$ 370.44

PH: 800-886-8675
 FX: 847-622-0454
 www.argos-tech.com

Distributed by: