

ARGOS INTRODUCES NEW INNOVATIONS IN CRYOWARE LABELING...

MICRYO DOTS AND STRIPS

A new range of self-adhesive labels designed specifically for low temperature use. This dedicated range of labels can be used for most applications to temperatures as low as -196 °C. They have a strong permanent adhesive (apply at +5 °C or more).

No need for separate labeling solutions for both your cryopreservation units and your ultra low freezers. They can even be used in incubators.

MICRYO DOTS AND STRIPS..... THE TRULY VERSATILE LABEL FROM ARGOS

The wide range of sizes allows you to label your whole experiment permanently. From microcentrifuge tube to cryovial, microtiter plate to storage box.

Micryo labels can be used on polypropylene, glass, polycarbonate and most metals.


ORDERING INFORMATION TUBES MICRYO DOTS AND STRIPS - ON ROLLS

Reference	Color	Dimensions	Quantity
MICRYO STRIPS			
L1271	White	33 x 13mm	roll of 1000 strips
L1286	Blue	33 x 13mm	roll of 1000 strips
L1287	Red	33 x 13mm	roll of 1000 strips
L1288	Green	33 x 13mm	roll of 1000 strips
L1289	Yellow	33 x 13mm	roll of 1000 strips
L1290	Purple	33 x 13mm	roll of 1000 strips

L1273	White	20 x 6mm	roll of 1000 strips
MICRYO DOTS			
L1274	White	9.5mm diameter	roll of 1000
L1275	White	9.5mm diameter	roll of 5000
L1292	5 colors	9.5mm diameter	roll of 5000*
*1000 of each Blue, Red, Green, Yellow, Purple			

L1276	White	19mm diameter	roll of 1000
L1277	White	19mm diameter	roll of 5000

ACTUAL SIZE AND COLORS





- No need to self-laminate
- Sleeker labels mean your tubes will still fit most storage systems (including cryocanes) without the label sticking or being torn off
- Easy to see, the bold range of 5 colors means that you can identify your samples while still in the box
- Unique adhesive formulation for that extra staying power
- Guaranteed not to fall off under normal use (providing user instructions are followed)
- A cost effective solution for your low temperature labeling needs





TUBEES® MICRYO LASER DOTS AND STRIPS

Designed for inkjet and laser printers, as well as other conventional devices. Supplied on sheets.



TUBEES MICRYO DOTS AND STRIPS - ON SHEETS

Reference	Color	Dimensions	Quantity
33 X 13MM, FOR USE ON 1.5/2ML CRYOVIALS, MICROCENTRIFUGE TUBES AND GLASS VIALS			
L1279	White	85 labels per sheet	pack of 20 sheets
L1280	Rainbow	85 labels per sheet	pack of 25 sheets*

24 X 13MM, FOR USE ON 0.5/1.5ML CRYOVIALS, PCR TUBES, MICROCENTRIFUGE TUBES AND GLASS VIALS

L1281	White	119 labels per sheet	pack of 20 sheets
L1282	Rainbow	119 labels per sheet	pack of 25 sheets**

67 X 25MM, FOR USE ON IDENTIFICATION RACKS, CRYOBOXES AND STORAGE BOXES

L1283	White	30 labels per sheet	pack of 20 sheets
-------	-------	---------------------	-------------------

38 X 6MM, FOR USE ON MICROTITER PLATES, PETRI DISHES AND STORAGE BOXES


L1284	White	156 labels per sheet	pack of 20 sheets
L1285	Rainbow	156 labels per sheet	pack of 25 sheets**

MICRYO LASER DOTS

SUITABLE FOR USE ON CRYOTUBES, GLASS VIALS AND MICROCENTRIFUGE TUBES

L1278	White, 19mm diameter	108 dots per sheet	pack of 20 sheets
L1291	White, 13mm diameter	192 dots per sheet	pack of 20 sheets

**15 sheets of each Blue, Red, Green, Yellow, Purple



Can also be used on cryoboxes



1141 East Main Street S104
 East Dundee, IL 60118
 Phone: 800-886-8675
 Fax: 847-783-0380
 www.argos-tech.com

