


# Are you using the right tools in your lab?



**It's okay to get caught short!**

# Introducing the *NEW* Argos Mini Serological Pipettes



- Measuring only 9" the Mini Pipet saves you 4" off standard pipettes
- Its petite size makes it ideal when aspirating in tissue culture hoods and tight spaces

- Polystyrene, gamma sterilized & pyrogen free
- Easy to read graduations from + 1 - 2%
- Unique filter plug stays put & eliminates lint fibers
- Color coded rings for easy identification



Fiber free paper wrap ensures quick peel apart opening

## Available in 3 *Mini* Sizes!

**PM5** - 5mL - 200/CS

**PM10** - 10mL - 150/CS

**PM25** - 25mL - 100/CS

\*Standard size pipettes also available

**ARGOS**  
TECHNOLOGIES

800-886-8675 / 847-622-0456  
[www.argos-tech.com](http://www.argos-tech.com)

*\*Available through your preferred distributor*

