




Pellet Mixer

Argos Technologies Pestle Grinder is a lightweight, battery powered mixer that can be used with disposable pestles and microtubes. Designed to homogenize tissue, cells or re-suspend protein and DNA pellets. Cell, tissue and pellets are completely homogenized within seconds by means of a high speed vortexing action. Battery powered unit greatly reduces hand fatigue.

Argos Pestle Grinder comes complete with 2AA batteries and pestle / mixer adapter.



Pestle & Tubes

Pestles are manufactured from polybutylene terephthalate (PBTP) in an ISO class 7 (FED-STD-209E class 10,000/M5.5) cleanroom. Both tubes and pestles are DNase, RNase and pyrogen free. The pestles feature a grooved grip handle that makes it easy to turn, even when wearing gloves.




Cat. No.	Description	Qty.
P7339-901	1.5mL Pestle	100/Pk.
P7338-901	0.5mL Pestle	100/Pk.
P7811-904	1.5mL Microtube	100/Pk.
P8185-904	0.5mL Microtube	100/Pk.
P9950-901	1.5mL Microtube & Pestle	100/Pk.
P9951-901	0.5mL Microtube & Pestle	100/Pk.
A0001	Pestle Motor Mixer	1/Ea.
A0002	Motor Adapter Replacement	1/Ea.

L-Shaped Cell Spreaders

Argos Technologies disposable, blue L-shaped, cell spreaders are ideal for spreading bacterial cultures evenly; the L-shape allows application on the entire petri dish or plate surface. Smooth spreading surface with slight upward return reduces the risk of damaging agar sample. Spreaders are manufactured of polypropylene and pre-sterilized, eliminating the need for flaming or autoclaving. Packaged in bags of 10 or individually in a peel back pouch.

Cat. No.	Description	Qty.
C4050	Cell Spreaders, PP, 10 per pack	500/Cs.
C4055	Cell Spreaders, PP Individually Wrapped Sterile	100/Cs.



see other side for more products!



Ph: 800-886-8675 / 847-622-0456
www.argos-tech.com