



polar safe


Reinforced lock and seal freezer storage boxes

Argos PolarSafe storage boxes are constructed of translucent polypropylene and can accommodate .5, 1.5 and 2.0mL microtubes as well as most brands of cryogenic vials. Argos brand boxes are slightly taller than most standard freezer boxes and can store tubes up to 5.3cm (2 1/8") in height. Boxes can withstand temperatures from -100°C to 121°C and will not degrade after repeated thaws like traditional cardboard boxes.


The reinforced seal, triple hinges and locking clasp keeps samples secure while the box is closed. Each box features an alpha-numeric grid on the exterior and interior, stacking rings and a matte surface suitable for printing and labeling.

Available with natural or fluorescent colored lids with natural bottoms for easy visual inspection. Black box is suitable for photo sensitive samples.

Dimensions: 14 X 6 X 14 cm. (5 1/2 X 2 3/8 X 5 1/2 in.)




Cat. No.	Description	Qty.
R3121	100-place box, black	5/pk.
R3125	100-place box, natural	5/pk.
R3126	100-place box, blue	5/pk.
R3127	100-place box, green	5/pk.
R2128	100-place box, orange	5/pk.
R3129	100-place box, purple	5/pk.
R3130	100-place boxes, mixed pack (one each of natural, blue, green, orange & purple)	5/pk.



Reinforced triple hinges prevent cracks or breaks, even when dropped



Bold, black alpha-numeric printing allows for accurate tracking of samples



Click fit lock ensures box will not open accidentally - large clasp for easy opening



Ph: 800-886-8675 • Fx: 847-622-0454
www.argos-tech.com

Distributed by: |