



Argos RotoFlex™

Benchtop Tube Rotator designed for a variety of applications!

The RotoFlex will accommodate a variety of tubes from 1.5 - 50mL by simply changing the rotisserie. No tools required! The mixing angle can be adjusted from 0-90° and each of the 4 tube panels can be aligned independently. The 360° rotation ensures thorough mixing of all samples at the same time while rotating at 22 rpm.

To change the rotation direction - simply touch the rotisserie paddle while in motion and the RotoFlex automatically begins to rotate in the opposite direction.

RotoFlex can be used in incubators & cold room settings. Unit is supplied with base unit and rotisseries (1 pack of each size).

**New
Item!**



Cat. No.	Description	Qty.
R2000	RotoFlex 120V	1/EA
R2001	RotoFlex 220V	1/EA
R2002	Rotisserie for .5 - .8mL tubes	2/PK
R2003	Rotisserie for 1.5 - 2mL tubes	2/PK
R2004	Rotisserie for 5 - 7mL tubes	2/PK
R2005	Rotisserie for 10 - 15mL tubes	2/PK
R2006	Rotisserie for 50mL tubes	1/PK

CE Approved / RoHS Compliant

Ideal for:

- Hybridization
- DNA Purification
- Sample mixing / resuspension
- Protein Binding
- Immunoprecipitation



Removable Inserts - easily snap in and out of place



Change inserts to accommodate varying tube sizes



Simply touch to engage auto reverse feature

| Distributed by: |

ARGOS

TECHNOLOGIES

Ph: 800-886-8675 • Fx: 847-622-0454
www.argos-tech.com

REV0511 LIT. NO. 9950