


15 & 50mL Polypropylene Storage Boxes

PolarSafe polypropylene freezer boxes can accommodate 15&50mL centrifuge tubes for long term storage; an excellent alternative to cardboard storage boxes. Boxes can be stored to -90°C (130°F) and are autoclavable to 122°C (251.6°F). (Dimensions: 5 3/4L x 5 3/4W x 4 7/8"H)



- Clear lids for easy visual inspection
- Easy-to-read numerical grid for tube reference and identification
- White top portion allows for easy labeling or marking
- Accommodates most 4" stainless steel freezer racks

| Catalog No. | Description | Qty. |
|-------------|--|--------|
| R0015 | 36-Place PP Freezer Box for 15mL Tubes | 2/Pack |
| R0050 | 16-Place PP Freezer Box for 50mL Tubes | 2/Pack |

15 & 50mL Cardboard Storage Box

Cardboard freezer boxes are constructed of sturdy, durable, white cardboard with a moisture repellant coating that protects the box from the affects of liquid nitrogen or condensation from refrigerators or freezers. Boxes feature small drain holes on the bottom for safe liquid nitrogen drainage or to allow excess moisture to escape.

Boxes and dividers sold separately.



Tubes sold separately


- Moisture repellant coating prevents water or liquid nitrogen from saturating box
- Easy to write on surface for identification
- Drain holes on bottom - allows for safe liquid nitrogen drainage and prevents excess moisture build-up

| Catalog No. | Description | Qty. |
|-------------|---|---------|
| R3014A | Cardboard Freezer Box (5 3/4 x 5 3/4 x 4 7/8 in.) | 1/Each |
| R3022 | 36-Place Divider for 15mL Tubes | 12/Pack |
| R3021 | 16-Place Divider for 50mL Tubes | 12/Pack |



15 & 50mL Polypropylene Tubes

Disposable polypropylene tubes and caps with recessed lid and sealing ring to prevent leaking. Frosted writing surface allows for permanent identification and labeling. Autoclavable to 121°C and can be stored in freezers to -80°C. Manufactured in ISO 9001:2000 and 100,000 grade cleanrooms. Sterile.



- Easy-to-read black graduations are accurate to $\pm 2\%$, 1mL increments (15mL) or 2.5mL increments (50mL)


- Available in cardboard racks or in resealable zip closure bags

- Can be spun up to 12,000g (15mL) & 15,000g (50mL)

- Conical Bottom

- Sterilized with gamma irradiation, RNase & DNase Free, DNA free, Non-pyrogenic

- Chloroform resistant graduations and writing area



Environmentally friendly, fully recyclable cardboard racks

| Catalog No. | Capacity | Packaging | Qty. |
|-------------|----------|----------------|-------------------|
| TC1500 | 15mL | Resealable Bag | 25/Bag, 500/Case |
| TC1500R | 15mL | Cardboard Rack | 25/Rack, 500/Case |
| TC5000 | 50mL | Resealable Bag | 25/Bag, 500/Case |
| TC5000R | 50mL | Cardboard Rack | 25/Rack, 500/Case |



Ph: 800-886-8675 • Fx: 847-622-0454
www.argos-tech.com

Distributed by: