



## Autoclave Deodorizers



These sterilizer deodorizers chemically absorb odor causing molecules, removing unwanted odors from your facility rather than masking them. Simply place one capsule in your autoclave or incubator with each load. The heat and moisture dissolves the capsule and releases a clean fresh scent.



Available in Apple, Heather and Lemon.


Cat. No.	Description	Qty
D6006C	Apple Deodorizer (40/Bottle)	6/CS
D6008C	Heather Deodorizer (40/Bottle)	6/CS
D6003C	Lemon Deodorizer (40/Bottle)	6/CS
D6005C	Apple Deodorizer (100/Bottle)	6/CS
D6007C	Heather Deodorizer (100/Bottle)	6/CS
D6004C	Lemon Deodorizer (100/Bottle)	6/CS

## Disposable Loops & Needles

High-impact loops and needles are smooth and flexible to facilitate uniform and smooth streaking without damaging the gel surface. Needles are straight and suitable for removal of single colony specimens. Packed sterile in safe, tamperproof, zip seal resealable bags. Disposable inoculating loops and needles do not need flaming, eliminating the risk of infection due to aerosol formation of pathogenic substances, as well as cross-contamination due to improper sterilization. Each case is supplied with a quality certificate confirming sterility and accuracy (Evans Blue Dye Method).

Color-coded for easy size identification.

Total length: 22.7cm.




Cat. No.	Description	Qty
<b>Inoculating Loop</b>		
L0001	1 Loop Blue, Sterile (25/Bag)	1000/CS
L0001-25	1 Loop Blue, Sterile (10/Bag)	1000/CS
L0010	10 Loop, Yellow Sterile (25/Bag)	1000/CS
L0010-25	10 Loop, Yellow Sterile (10/Bag)	1000/CS
<b>Inoculating Needle</b>		
L0000	Needle, White Sterile (25/Bag)	1000/CS
L0000-25	Needle, White Sterile (10/Bag)	1000/CS

**Argos Technologies, Inc.**

1551 So. Scottsdale Ct. S200  
Elgin, IL 60123

Ph: 800-886-8675  
FX: 847-622-0454  
[www.argos-tech.com](http://www.argos-tech.com)



## Inoculating Loops & Needles:

Argos brand Inoculating Loops are made from high quality twisted nichrome wire. Loops have been designed to slide along agar surfaces without cutting or tearing. Argos Inoculating Loops cool rapidly and are less expensive than platinum. Loops are supplied without handles and are available calibrated or non-calibrated. Calibrated Loops are supplied with calibration certificates based on ASM/APHA Evans Blue Dye test protocols.

## Inoculating Loop Holders:

Argos Inoculating Loop Holders are made from nickel plated brass. Handles are insulated to protect the hand during use. Loop Holders are available in two sizes: 6" and 8".

## Inoculating Handles and Loops:

Reusable Inoculating Loops with handles are used to select, inoculate or transfer microbiological samples. Loops are made of high quality twisted nichrome wire and attached to a light-weight aluminum handle. Handles are insulated to protect the hand during use. Loops have been designed to slide along agar surfaces without cutting or tearing. Loops and Handles cool rapidly after sterilization and are offered calibrated (1µl & 10µl) or non-calibrated. Calibrated loops are supplied with calibration certificates based on ASM/APHA Evans Blue Dye test protocols. Non-calibrated Loops and Handles are available as small, medium or large and are color coded for easy identification.

## Inoculating Needles:

Reusable Inoculating Needles with handles to select, inoculate or transfer microbiological samples. Needles are made from high quality nichrome wire that cools rapidly after sterilization. Attached to a light-weight aluminum handle, insulated to protect the hand during use. Teasing Needles are offered in 6" straight (white) or 6" angled (blue).

Cat. No.	Description	Qty
<b>Loops without Handles</b>		
L1102	Small Loop / 2.26mm ID	25/PK
L1103	Medium Loop / 2.91mm ID	25/PK
L1104	Large Loop / 3.57mm ID	25/PK
<b>Calibrated Loops without Handles</b>		
L1100	Calibrated Loop, 10ul / 4mm ID	5/PK
L1101	Calibrated Loop, 1ul / 1.45mm ID	5/PK
<b>Streaker Loops with Handles</b>		
L1120	Small Streaker, Blue / 2.26mm ID	5/PK
L1121	Medium Streaker, Green / 2.91mm ID	5/PK
L1122	Large Streaker, Red / 3.57mm ID	5/PK
<b>Calibrated Streaker Loops with Handles</b>		
L1123	Calibrated Streaker 1ul / 1.45mm ID	5/PK
L1124	Calibrated Streaker 10ul / 4mm ID	5/PK
<b>Needles</b>		
L1125	Micro , Clear	5/PK
<b>Loop Holders</b>		
L1110	Brass Loop Holder, 6"	1/EA
L1111	Brass Loop Holder, 8"	1/EA