

# Expand Your Limits!

The "just right" size tube

## 5mL Graduated Prep Tube



The Argos 5mL tube is designed for samples too large for standard microcentrifuge tubes (1.5 - 2.0mL) and too small to pipette without contamination from a 15mL centrifuge tube. The easy-to-open cap allows for single handed operation without the risk of hand fatigue while still maintaining sample integrity.

Produced from medical grade polypropylene co-polymer, all tubes provide clear viewing of samples. Tubes are produced in a class 100,000 clean room and are RNase, DNase, pyrogen free. Tubes will withstand autoclaving (121° C for 20 minutes at 15 psi), suitable for freezing to -90° C (130° F) and can be safely centrifuged up to 14,000 rpm. A frosted marking area on the tube cap is provided for sample tagging and identification.

- Available in Natural, Blue, Rose, Green, Yellow and photo-resistant Black, sold in packs of 200.

- Gamma irradiated sterile tubes are available in Natural, sold in packs of 200 (20/pack, 10/packs per case)

### Ordering Information:

T2076A	Natural Non-sterile / 200 per pack	T2079A	Green Non-sterile / 200 per pack
T2077A	Blue Non-sterile / 200 per pack	T2080A	Rose Non-sterile / 200 per pack
T2078A	Yellow Non-sterile / 200 per pack	T2081A	Black Non-sterile / 200 per pack
		T2076S-CA	Natural - Sterile / 20 per pack, 10 pack/case

## FlexiFuge Centrifuge

FlexiFuge is a versatile, personal micro-centrifuge for quick, low speed centrifugation. Whisper quiet motor offers two speeds (2000g or 5000g / 6,000 or 9,200 RPM ) depending on application. FlexiFuge can be operated in continuous mode, pulsed for a short burst of rapid centrifugation, or timed to spin from 1-15 minutes. Regardless of setting, rotors stop once the lid release button is pushed. A ring of LED lights illuminate when the unit is active and can be viewed through the transparent lid. Color can be changed (red, green, blue) or turned off by a switch under the unit.

### Ordering Information:

C1000	FlexiFuge Centrifuge 120V (includes 1 of each rotor) / 1 each
C0300	1.5mL Rotor / 1 each
C0301	PCR Strip Rotor / 1 each
C0302	5mL Rotor / 1 each
C0303	Adapters for .5mL Tubes, 8 each



FlexiFuge is supplied with an 8-Place fixed angle rotor for 1.5mL tubes, 4-Place rotor for 5mL tubes, rotor for (4), 8 x.2mL PCR tube strips and 8 adapters for .5mL tubes.