



# MiniVac

## Simple, safe, cost effective waste removal system

The Argos MiniVac waste collection system was designed with safety, convenience and bench-top space in mind. The MiniVac system replaces the need for glass waste containers; eliminating the risk of breakage and contamination. MiniVac can be used as a stand-alone waste collection system or as an intermediate overflow vessel. The unit consists of a disposable 1L polystyrene flask, plastic reservoir basin, polypropylene cap with barbed fittings and tubing. One tubing set is connected to primary or secondary source and one is used to aspirate the liquid. MiniVac flask and cap sit in a clear PETG plastic reservoir basin that collects liquid should flask crack or leak. Full flasks can be transferred safely by using supplied transfer cap. Small footprint (7.5 in.); easy clean up with 70% isopropyl alcohol.



**New  
Safety Item!**

- Safely aspirate & transfer waste
- Clear safety reservoir collects liquid in case flask cracks or leaks
- Use with in-house vacuum
- Polypropylene barbed cap fittings eliminate risk of broken glass rod fittings
- Use as stand-alone waste collection system or overflow vessel

Cat. No.	Description	Qty.
MV000	MiniVac Waste Aspiration System, Complete	1/Set
MV001	Polypropylene barbed cap	1/EA
B0001	Polystyrene bottle with transfer cap	1/EA
B0000	Reservoir basin, PETG plastic	1/EA
EV331	Tubing 700mm	1/EA

| Distributed by: |

# ARGOS

TECHNOLOGIES

Ph: 800-886-8675 • Fx: 847-622-0454  
www.argos-tech.com

LIT. NO. 9951L-A / REV0212